Quiz 1 BT 601 Time: 45 min Total Marks: 10

[Show all steps and No part marking for numerical problems]

- when the retention time is 10.50 min and the half-base width is 1.2 min and 1. Calculate the number of theoretical pallets (N) and the plate height (H), the column length is 25 cm.
 - 2. a) Show the relation (expression) of Kn and Kp.
- gm of the b) Calculate the KP of a compound, when 10 gm of the compound present in 300 ml water was shaken with 20 ml of acetone and 8 compound was transferred to acetone layer only.
 - 3. Write the distribution of 1000 molecules of a compound in 5 theoretical plates with K_p=0.1 [Consider whole number for distribution].
- 4. a) Write the expression related to Kp and Ve based on the column parameters
- and "Include" b) Modify the above expression for totally "exclude" conditions.
 - 5. a) Show the chromatogram of different Mw proteins in a gel filtration chromatography?
 - b) Show the plot to determine unknown Mw of a protein.

[Marks: 2+2+2+2= 10]

